



PO Box 788
Belmont, NC 28012
+1.800.732.8484 F: +1.888.843.4384
Reagents.com

info@Reagents.com
Phone: 704.554.7474
Toll Free: 800.732.8484
Fax: 888.843.4384

Specification Sheet

Product: Lead Carbonate, ACS Reagent

Catalog Number: C1227000

Shelf Life: 5 Years

CAS Number: 598-63-0

Recommended Storage: 15 to 25 °C (59 to 77 °F)

Test	Specification
Nitrate and Nitrite (as NO ₃)	To pass test
Sodium (Na)	0 .05 % MAXIMUM
Chloride (Cl)	0 .002 % MAXIMUM
Potassium (K)	0 .02 % MAXIMUM
Calcium (Ca)	0 .01 % MAXIMUM
Iron (Fe)	0 .005 % MAXIMUM
Zinc (Zn)	0 .003 % MAXIMUM
Cadmium (Cd)	0 .002 % MAXIMUM
Insoluble in Dilute Acetic Acid	0 .02 % MAXIMUM