

## FKBP3 Polyclonal Antibody

Catalog number: 11755-1-AP


Size: 57 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

FKBP3; FKBP-25; PPIase



HeLa cells were subjected to SDS PAGE followed by western blot with 11755-1-AP(FKBP3 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells, A431 cells, human brain tissue, PC-3 cells
Calculated FKBP3 MW:	224aa,25kd
Observed FKBP3 MW:	25kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag2351
GenBank accession number:	BC020809
Gene ID (NCBI):	2287
Full name:	FK506 binding protein 3, 25kDa

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.