

A4GALT Polyclonal Antibody

Catalog number: 12392-1-AP

Size: 30 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

A4GALT; A4GALT, A4GALT1,

CD77 synthase, Gb3 synthase



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 12392-1-AP(A4GALT Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated A4GALT MW:	353aa, 41 kDa
Observed A4GALT MW:	50-52 kDa
Positive WB detected in	RAW 264.7 cells, K-562 cells, Raji cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3069
GenBank accession number:	BC017068
Gene ID (NCBI):	53947
Full name:	Alpha 1,4-galactosyltransferase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.