

SALL2 Polyclonal Antibody

Catalog number: 12679-1-AP

Size: 39 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

SALL2; FLJ10414; FLJ55746;

HSAL2; KIAA0360; ZNF795;

p150(Sal2)



HEK-293 cells were subjected to SDS PAGE followed by western blot with 12679-1-AP(SALL2 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HEK-293 cells
Calculated SALL2 MW:	1007aa,105kd
Observed SALL2 MW:	105kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3367
GenBank accession number:	BC024245
Gene ID (NCBI):	6297
Full name:	Sal-like 2 (Drosophila)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.