

## FAM19A5 Polyclonal Antibody

Catalog number: 13948-1-AP


Size: 57 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

FAM19A5; Chemokine like protein TAFA 5, FAM19A5, Protein FAM19A5, QLLK5208, TAFA 5, TAFA5, UNQ5208



Immunohistochemistry of paraffin-embedded human brain tissue slide using 13948-1-AP (FAM19A5 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using 13948-1-AP (FAM19A5 Antibody) at dilution of 1:50 (under 40x lens)

### Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated FAM19A5 MW:	14 kDa
Observed FAM19A5 MW:	
Positive IHC detected in	Human brain tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5058
GenBank accession number:	BC039396
Gene ID (NCBI):	25817
Full name:	Family with sequence similarity 19 (chemokine (C-C motif)-like), member A5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.