Constant Temperature Water Bath


BS Series

Operating temperature range RT+5°C~Water boiling point

Temperature distribution accuration

±5°C (at 70°C) Bath capacity

Standard equipped with overheat prevention device, electric leakage breaker.



BS200

4.7L BS200

91

BS400

 4 types of accessorial baskets. Maximum 5L beaker can be set after removing the lid

14L

BS660

11L BS600

BS400

Comes with 4 container mounting clamps and cooling water distributor (optional item)

BS600

 Up to 6pcs. of containers can be used simultaneously

BS660

Fits small and large containers. Up to 6pcs. can be used simultaneously. Sink level can be adjusted

Control Panel



For different flask size

Specifications

Features

Easy to change water level by adjusting over flow

C14122

Removable lid



Sturdy die-cast clamps

BS400 Model **BS200 BS600 BS660** Operating temperature range RT+5°C~Water boiling point *1 Temperature adjustment accuracy ±3°C (at 70°C) Temperature distribution accuracy ±5°C (at 70°C) Max. temperature reaching time ~35min. ~60min. ~50min. Stainless steel Internal tank material Heater Copper Pipe Heater (Nickel Plated) 0.9kW 1.4kW 1.3kW Safety device Overheat prevention system (Liquid expansion type, Self sustaining relay, Alarm lamp), Electric leakage breaker with overcurrent protection L790 x W150 x D102mm Internal dimension ø200 x D149mm L328 x W302 x D99mm L502 x W302 x D99mm Effective internal dimension ø186 x D115mm L298 x W260 x D64mm L776 x W140 x D73mm L468 x W260 x D66mm External dimension *3 L390 x W438 x H214mm L864 x W288 x H215mm L300 x W364 x H216mm L565 x W437 x H214mm Bath capacity ~4.7L ~9L ~11L ~14L Power source (50/60Hz) AC115V, 11.5A AC115V, 8A AC115V, 12.5A AC220V, 4.5A AC220V, 6.5A AC220V, 6A *2 Opening size, quantity ø115mm, 4pcs. ø115mm, 6pcs ~16kg Weight ~5kg ~9.5kg ~12ka Accessories Ring (excluding BS200), lid, clamp, stand, clamp bar (excluding BS200), connector (for stand), thermometer (0~+100°C alcohol type), bottom plate and basket (for BS200)

*1. Ambient temperature might impact the boiling point. Conditions: RT23°C±5°C, humidity 65%±20% (no load).
*2. 4 lids. Suitable for erlenmeyer flask up to ø96 x D50mm 300ml, ø56 x D40mm 50ml, ø40 x 30mm.

*3. Protrusions not included. Power cord lenght 2m.

Optional items

Product code	Description	Suitable models
222081	Cooling water distributor (supply/drain tube ø12mm, caliber of spout ø8.5mm)	BS400
222086	Cooling water distributor (supply/drain tube ø16mm, caliber of spout ø8.5mm)	BS600

Includes distribution pipe 2pcs., connector 4pcs., 2pcs.=1set (supply/drain). One connected to water faucet will inject water to 4 cooling water distributors.



Cooling water distributor for BS400