

CYTSB Polyclonal Antibody

Catalog number: 13851-1-AP

Size: 48 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

CYTSB; FLJ36955; HCMOGT-1;

NSP; SPECC1



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 13851-1-AP(CYTSB Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse kidney tissue, HeLa cells
Calculated CYTSB MW:	119kd
Observed CYTSB MW:	119kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4822
GenBank accession number:	BC050058
Gene ID (NCBI):	92521
Full name:	Cytospin B

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.