




### Specifications

|                                  |  |
|----------------------------------|--|
| Electrical requirements          | 208-240VAC 50/60Hz single phase, 2500 watts  |
| Chamber dimensions               | 22 x 16 x 22" (WxDxH)                        |
| Unit dimensions                  | 30 x 26 x 34" (WxDxH)                        |
| Shipping dimensions              | 33 x 32 x 41" (WxDxH)                        |
| Temperature range                | Ambient temperature to 300°C                 |
| Temperature uniformity           | +/- 2°F                                      |
| Temperature controller precision | 0.1°F  |
| Timer range                      | 1 to 9999 minutes                            |
| Racks included                   | 2 (additional racks available at extra cost) |
| Maximum rack capacity            | 9 (height between each rack is 2 inches)     |
| Exhaust port                     | 2" diameter x 1, adjustable                  |
| Shipping weight                  | 275 Lb                                       |
| Access port                      | 2" diameter with cover, one                  |
| Certification                    | CE, UL/CSA (optional)                        |

### Standard package

| Part description                  | Quantity | Part image  |
|-----------------------------------|----------|---|
| AI 4.8 cu ft forced air oven      | 1 pc     |  |
| Stainless steel rack and brackets | 2 sets   |  |
| User's manual                     | 1 pc     |  |