

Whatever Your Labeling Need, We Have the Perfect Pen or Tool for You!



	A	B	C	D	E	F
	Wet Surface Pen	Tech Pen	Belpen™	Autoradiology Pens	Captor™ Pen Holder	Engraving Tool
Unique Feature	Writes on most surfaces	Oil based marking pen resists washing, weathering, and mild acids and alkalis	Quick drying waterproof markers	Quick, accurate X-ray film marking	Cap and uncup pens with one hand	Portable and permanently engraves in seconds
Marker			X	X		
Inscriber						X
Mark Size	Fine Point	Medium Point	Medium Point	Fine Point	N/A	Medium Point
Temperature Resistant	-34°C to 121°C (-30°F to 250°F)	-46°C to 65°C (-50°F to 150°F)				
Withstands Autoclaving		X	X			
Waterproof	X	X	X			
Smudge Proof			X			
Available in Colors		X	X			
Pocket Clip	X		X	X		
Operates on Two 'AA' Batteries						X (not included)
WRITES ON:						
Ceramics		X				X
Cloth						
Glassware		X	X			X
Greasy Surfaces	X					
Hard to Mark Surfaces	X					
Label Tape	X	X				
Metal		X				X
Paper	X	X				
Plastics		X	X			X
Radiographical Films				X		
Stone / Bone		X	X			X
Teflon® Fluoropolymer Resins						X
Wet Surfaces	X					
Wood		X				X

Whatever Your Labeling Need, We Have the Perfect Pen or Tool for You!

Complete collection of pens and tools that work on virtually all materials:

- Wet Surface Pen writes clearly on wet, greasy and other hard to mark surfaces. Use it under water, upside down and in the harshest of temperatures
- Tech Pens write on glass, metal, ceramics, plastics, paper
- Belpen™ Markers contain permanent and quick drying ink that is waterproof and smudge proof
- Autoradiography Pens are excellent for marking substrates that are to be radiographed
- Captor™ Pen Holder keeps a pen handy and easy to uncap with only one hand
- Battery Powered Engraving Tool permanently engraves in seconds

Use the chart on the reverse side to see which marking tools are best suited to your surface.

Product	Description	Unit of Measure	Thomas No.
Wet Surface Pen	Writes on Wet Surfaces, under Water and at any Angle, Black Ink	Each	1226C14
Tech Pens	Write on Most Surfaces, Green	Each	1222X48
Tech Pens	Write on Most Surfaces, Black	Each	1222X43
Tech Pens	Write on Most Surfaces, White	Each	1222X45
Tech Pens	Write on Most Surfaces, Orange	Each	1222X49
Tech Pens	Write on Most Surfaces, Blue	Each	1222X47
Tech Pens	Write on Most Surfaces, Yellow	Each	1222X46
Tech Pens	Write on Most Surfaces, Red	Each	1222X44
Belpen™ Markers	Quick Drying Waterproof Ink, Black	3 per pack	1226C12
Belpen™ Markers	Quick Drying Waterproof Ink, Blue	3 per pack	1226C13
Autoradiography Pens with Handy Clip	Quick, Accurate Film Marking Normal Energy Level, Non-Radioactive Marker	Each	1226C06
Autoradiography Pens with Handy Clip	Quick, Accurate Film Marking High Energy Level, Non-Radioactive Marker	Each	1226C07
Captor™ Pen Holder	Cap and Uncap Pens with One Hand	3 per bag	1222X42
Battery Powered Engraving Tool	Permanently Engraves in Seconds, on Many Surfaces Including Glass, Plastic and Metal	Each	1221Z88
Engraving Tool Stand	Holds Engraving Tool When Not In Use	Each	1221Z89

