



PO Box 788
Belmont, NC 28012
+1.800.732.8484 F: +1.888.843.4384
Reagents.com

Specification Sheet

Product: HPLC Water
Also meets ACS specifications

Catalog Number: C4346000

Shelf Life: 2 Years

CAS Number: 7732-18-5

Recommended Storage: 15 to 25 °C (59 to 77 °F)

Test	Specification
Appearance	Clear, colorless
Color (APHA)	5 MAXIMUM
Conductivity at 25°C	1.5 µS/cm MAXIMUM
LC Suitability (Absorbance)	To pass test
Residue after Evaporation	1 ppm MAXIMUM
Substances Reducing MnO ₄	To pass test
Fluorescence at 450nm	1 ppb MAXIMUM
Absorbance at 200 nm	0 .01 AU/cm MAXIMUM
Absorbance at 210 nm	0 .01 AU/cm MAXIMUM
Absorbance at 250 nm	0 .005 AU/cm MAXIMUM
Absorbance at 400 nm	0 .005 AU/cm MAXIMUM