

Product Name:	ABS VS 7
Description:	<p>ABS VS Class Vapor Shipper Total Vial Capacity (6/Cane): 500 Number of Racks: 1 Boxes per Rack: 5 Total Bag Capacity (250 mL): 10 Liquid Nitrogen Capacity: 8.5 Liters Static Evaporation Rate (Liters/day): 0.85 Static Holding Time (days): 10 Neck Diameter (mm): 215.9 Exterior Height (mm): 546.1 Exterior Diameter (mm): 368.3 Optional protective shipping container ensure upright shipping High strength neck tube reduces liquid nitrogen loss Lid with locking provision Hydrophobic liquid nitrogen absorbent system 2-hour charging capability One-year parts and labor warranty One-year vacuum warranty</p>