








 BIOFIL ®	Document Number	Revision	Page	Document Number
	ATT003	A0	1 of 2	12-04-2020
PRODUCT DATA SHEET		Robotic Tips		



Dimensions (Unit: mm)



Robotic Tips

Description: Tecan-Style Robotic tips, Rack

Purpose: Accurate liquid handling

Materials

Tips: PP (Polypropylene)

Rack box: PP (Polypropylene)

Features

- JET Robotic Tips are made from pure Polypropylene
- Available in both black conductive tips and clear non-conductive tips
- Black conductive tips: liquid level sensing
- Compatibility: Tecan Genesis Freedom®, Freedom Evo® and Miniprep with LiHa® workstations
- Tip format: 96 tips/Rack
- Tip volume range selection: 20µL, 50µL, 200µL and 1000µL
- Available with or without aerosol-resistant filters
- Low Retention/Standard surface
- DNase/RNase-free and non-pyrogenic
- Endotoxin-Safe
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Robotic Tips

Cat. No.	Max Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
ATT000020	20μL	Natural	Standard	No	No	96 tips/rack, 24 racks/case.
AMT000020			Low Retention	No	No	
ATT001020			Standard	No	Yes	
AMT001020			Low Retention	No	Yes	
ATT101010	10μL		Standard	Yes	Yes	
AMT101010			Low Retention	Yes	Yes	

 BIOFIL ®	Document Number	Revision	Page	Document Number
	ATT003	A0	2 of 2	12-04-2020
PRODUCT DATA SHEET		Robotic Tips		

Robotic Tips

Cat. No.	Max Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
AUT000020	20μL	Black	Standard	No	No	96 tips/rack, 24 racks/case. (Conductive black)
ANT000020			Low Retention	No	No	
AUT001020			Standard	No	Yes	
ANT001020			Low Retention	No	Yes	
AUT101010	10μL		Standard	Yes	Yes	
ANT101010			Low Retention	Yes	Yes	
ATT000050	50μL	Natural	Standard	No	No	96 tips/rack, 24 racks/case.
AMT000050			Low Retention	No	No	
ATT001050			Standard	No	Yes	
AMT001050			Low Retention	No	Yes	
ATT101050			Standard	Yes	Yes	
AMT101050			Low Retention	Yes	Yes	
AUT000050		Black	Standard	No	No	96 tips/rack, 24 racks/case. (Conductive black)
ANT000050			Low Retention	No	No	
AUT001050			Standard	No	Yes	
ANT001050			Low Retention	No	Yes	
AUT101050			Standard	Yes	Yes	
ANT101050			Low Retention	Yes	Yes	
ATT000200	200μL	Natural	Standard	No	No	96 tips/rack, 24 racks/case.
AMT000200			Low Retention	No	No	
ATT001200			Standard	No	Yes	
AMT001200			Low Retention	No	Yes	
ATT101200			Standard	Yes	Yes	
AMT101200			Low Retention	Yes	Yes	
AUT000200		Black	Standard	No	No	96 tips/rack, 24 racks/case. (Conductive black)
ANT000200			Low Retention	No	No	
AUT001200			Standard	No	Yes	
ANT001200			Low Retention	No	Yes	
AUT101200			Standard	Yes	Yes	
ANT101200			Low Retention	Yes	Yes	
ATT000000	1000μL	Natural	Standard	No	No	96 tips/rack, 16 racks/case.
AMT000000			Low Retention	No	No	
ATT001000			Standard	No	Yes	
AMT001000			Low Retention	No	Yes	
ATT101000			Standard	Yes	Yes	
AMT101000			Low Retention	Yes	Yes	
AUT000000		Black	Standard	No	No	96 tips/rack, 16 racks/case. (Conductive black)
ANT000000			Low Retention	No	No	
AUT001000			Standard	No	Yes	
ANT001000			Low Retention	No	Yes	
AUT101000			Standard	Yes	Yes	
ANT101000			Low Retention	Yes	Yes	

Drawn	Approved
Date 04 Dec 2020	Date 04 Dec 2020
Name Alex Fang	Name James