

## CYB5R4 Polyclonal Antibody

Catalog number: 12626-1-AP

Size: 25 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

CYB5R4; NCB5OR; RP4-

676J13.1; dJ676J13.1



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 12626-1-AP (CYB5R4 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	BxPC-3 cells, HeLa cells
Calculated CYB5R4 MW:	55.8kd
Observed CYB5R4 MW:	55.8kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3318
GenBank accession number:	BC025380
Gene ID (NCBI):	51167
Full name:	Cytochrome b5 reductase 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.