

SUV39H2 Polyclonal Antibody

Catalog number: 11338-1-AP

Size: 31 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

SUV39H2; FLJ23414, H3 K9

HMTase 2, KMT1B, Lysine N

methyltransferase 1B, Su(var)3

9 homolog 2, SUV39H2



HEK-293 cells were subjected to SDS PAGE followed by western blot with 11338-1-AP(SUV39H2 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Cited species:	Human
Calculated SUV39H2 MW:	53 kDa
Observed SUV39H2 MW:	53 kDa
Positive WB detected in	HEK-293 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1882
GenBank accession number:	BC007754
Gene ID (NCBI):	79723
Full name:	Suppressor of variegation 3-9 homolog 2 (Drosophila)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.