

## SCARF1 Polyclonal Antibody

Catalog number: 13702-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

SCARF1; Acetyl LDL receptor,

KIAA0149, SCARF1, SREC



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13702-1-AP(SCARF1 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated SCARF1 MW:	830aa,87 kDa
Observed SCARF1 MW:	86 kDa
Positive WB detected in	Mouse lung tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4621
GenBank accession number:	BC039735
Gene ID (NCBI):	8578
Full name:	Scavenger receptor class F, member 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.