

V-PETTE

STARTER KIT



- SET OF 4 COLOR PIPETTORS AND ACCESSORIES
- VOLUME RANGE FROM 0.5 μ L UP TO 1000 μ L
- DURABLE, UV AND CHEMICALLY RESISTANT PIPETTORS
- UNIVERSAL 4-POSITION PIPETTE STAND
- UNIVERSAL PIPETTE TIP COMPATIBILITY

V-PETTE STARTER KIT ORDERING INFORMATION

CAT. No.	DESCRIPTION
P3946-SK4	SINGLE CHANNEL VARIABLE VOLUME PIPETTES: VP10 , VP20 , VP200 , VP1000 4-POSITION PIPETTE STAND 3 RACKS OF PIPETTE TIPS INSTRUCTION MANUAL ACCESSORIES

TECHNICAL INFORMATION

MODEL	VOLUME [μ L]	ACCURACY [%]	PRECISION [%]
VP10	0.5 - 10	MIN \pm 4.0 - MAX \pm 0.5	MIN \pm 2.8 - MAX \pm 0.4
VP20	2 - 20	MIN \pm 3.0 - MAX \pm 0.8	MIN \pm 1.5 - MAX \pm 0.3
VP200	20 - 200	MIN \pm 1.2 - MAX \pm 0.6	MIN \pm 0.6 - MAX \pm 0.2
VP1000	100 - 1000	MIN \pm 1.6 - MAX \pm 0.6	MIN \pm 0.4 - MAX \pm 0.15

 PZ HTL S.A., Daniszewska 4
03-230 Warsaw, POLAND

All rights reserved. Product described in this brochure is subject to availability and technical modification.

Errors excepted. PZ HTL S.A. reserves the right to improve, enhance or otherwise modify its products without prior notification. © 2014 PZ HTL S.A.