

DH5a-lysate Polyclonal Antibody

Catalog number: 10559-1-Ig

Size: 105 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DH5a lysate

Applications

Tested applications: ELISA, WB
Species specificity: Bacterial; other species not tested.
Calculated DH5a-lysate MW: 1-300 kDa
Observed DH5a-lysate MW:
Positive WB detected in Total bacteria lysates of DH5 alpha
Recommended dilution: WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Bacterial Lysates
GenBank accession number:
Gene ID (NCBI):
Full name: Undefined

Product information

Purification method: Protein A purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.