



Portable Freezer

Wide temperature range:

+10°C to -40°C

Stores up to

7,000 vaccine doses

IN STOCK Thomas No. 1213C23



**25 L (0.9 cu. ft.)
capacity**



ThomasSci.com 833.544.SHIP (7447)
CustomerService@ThomasSci.com

Connect With Us:



PRODUCT SPECIFICATIONS

Capacity	25 Litres (0.9 cu. ft.)
Temperature Range	+10°C to -40°C, 1°C increments
Power Source Voltage	DC 12V
Exterior Dimensions (W x D x H)	27.36" x 13.78" x 18.11"
Interior Dimensions (W x D x H)	13.19" x 8.86" x 13.39"
Noise	55 dB(A)
Power Consumption	48 VA
Weight	15.5 kg
Includes	AC adapter and vehicle adapter



Contact your
Representative to place
an order today!

Canadawide
SCIENTIFIC

Thomas
Scientific

ThomasSci.com 833.544.SHIP (7447)
CustomerService@ThomasSci.com

Connect With Us:

