

# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## 1. Identification

### 1.1. Product identifier

<b>Product Identity</b>	SDNA-1000
<b>Alternate Names</b>	SDNA-1000

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

<b>Intended use</b>	DNA/RNA preservation solution.
---------------------	--------------------------------

### 1.3. Details of the supplier of the safety data sheet

<b>Company Name</b>	Spectrum Solutions 12248 S Lone Peak PKWY Suite 106, Draper, UT 84020
---------------------	---

### Emergency

<b>CHEMTREC (USA)</b>	(800) 424-9300
<b>24 hour Emergency Telephone No.</b>	+1 650-209-8338
<b>Customer Service: Spectrum Solutions</b>	+1 650-209-8338

## 2. Hazard(s) identification

### 2.1. Classification of the substance or mixture

Flam. Liq. 3;H226	Flammable liquid and vapor.
Acute Tox. 4;H302	Harmful if swallowed.
Acute Tox. 5;H313	May be harmful in contact with skin. (Note: this GHS classification was not adopted by US OSHA).
Skin Irrit. 3;H316	Causes mild skin irritation. (Note: this GHS classification was not adopted by US OSHA).

### 2.2. Label elements



### Warning

H226 Flammable liquid and vapor.  
 H302 Harmful if swallowed.  
 H313 May be harmful in contact with skin.  
 H316 Causes mild skin irritation.

### [Prevention]:

P210 Keep away from heat, sparks, open flames, and other ignition sources - No smoking.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

- P233 Keep container tightly closed.  
 P240 Ground, bond container and receiving equipment.  
 P241 Use explosion-proof electrical, ventilating, light, equipment.  
 P242 Use only non-sparking tools.  
 P243 Take precautionary measures against static discharge.  
 P264 Wash thoroughly after handling.  
 P270 Do not eat, drink or smoke when using this product.  
 P280 Wear protective gloves, eye protection, face protection.

## [Response]:

- P301+312 IF SWALLOWED: Call a POISON CENTER, doctor or physician if you feel unwell.  
 P303+361+353 IF ON SKIN (or hair): Remove, take off immediately all contaminated clothing. Rinse skin with water, shower.  
 P312 Call a POISON CENTER, doctor or physician if you feel unwell.  
 P330 IF SWALLOWED: Rinse mouth.  
 P332+313 IF SKIN IRRITATION OCCURS: Get medical advice or attention.  
 P370+378 In case of fire: Use extinguishing media listed in section 5 of SDS for extinction.

## [Storage]:

- P403+235 Store in a well ventilated place. Keep cool.

## [Disposal]:

- P501 Dispose of contents or container in accordance with local and national regulations.

## 3. Composition/information on ingredients

This product contains the following substances that present a hazard within the meaning of the relevant State and Federal Hazardous Substances regulations.

Ingredient/Chemical Designations	Weight %	GHS Classification	Notes
Thiocyanic acid, compound with guanidine (1:1) CAS Number: 0000593-84-0	25 - 50	Acute Tox. 4;H302 Acute Tox. 4;H312 Acute Tox. 4;H332 Aquatic Chronic 3;H412	----
Ethanol CAS Number: 0000064-17-5	10 - 25	Flam. Liq. 2;H225	----
Tris (hydroxymethyl) aminomethane CAS Number: 0000077-86-1	1 - 5	Skin Irrit. 2;H315 Eye Irrit. 2;H319 STOT SE 3;H335	----

In accordance with paragraph (i) of §1910.1200, the specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

\*PBT/vPvB - PBT-substance or vPvB-substance.

The full texts of the phrases are shown in Section 16.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## Section 4. First aid measures

### 4.1. Description of first aid measures

<b>General</b>	In all cases of doubt, or when symptoms persist, seek medical attention. Never give anything by mouth to an unconscious person.
<b>Inhalation</b>	Remove to fresh air, keep patient warm and at rest. If breathing is irregular or stopped, give artificial respiration. If unconscious, place in the recovery position and obtain immediate medical attention. Give nothing by mouth.
<b>Eyes</b>	Irrigate copiously with clean water for at least 15 minutes, holding the eyelids apart and seek medical attention.
<b>Skin</b>	Remove contaminated clothing. Wash skin thoroughly with soap and water or use a recognized skin cleanser.
<b>Ingestion</b>	If swallowed obtain immediate medical attention. Keep at rest. Do NOT induce vomiting.

### 4.2. Most important symptoms and effects, both acute and delayed

<b>Overview</b>	Treat symptomatically. Exposure to solvent vapor concentrations from the component solvents in excess of the stated occupational exposure limits may result in adverse health effects such as mucous membrane and respiratory system irritation and adverse effects on the kidneys, liver and central nervous system. Symptoms include headache, nausea, dizziness, fatigue, muscular weakness, drowsiness and in extreme cases, loss of consciousness.  Repeated or prolonged contact with the preparation may cause removal of natural fat from the skin resulting in dryness, irritation and possible non-allergic contact dermatitis. Solvents may also be absorbed through the skin. Splashes of liquid in the eyes may cause irritation and soreness with possible reversible damage. Check section 2.2 (GHS Label Elements) for further details.
<b>Skin</b>	May be harmful in contact with skin. Causes mild skin irritation.
<b>Ingestion</b>	Harmful if swallowed.

## Section 5. Fire-fighting measures

### 5.1. Extinguishing media

Use extinguishing agents appropriate for surrounding fire: Regular dry chemical, carbon dioxide (CO<sub>2</sub>), regular foam

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition: Decomposition products may include carbon dioxide, carbon monoxide and nitrogen oxides.

Keep away from heat, sparks, open flames, and other ignition sources - No smoking.

Keep container tightly closed.

Ground, bond container and receiving equipment.

Use explosion-proof electrical, ventilating, light, equipment.

Use only non-sparking tools.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

Take precautionary measures against static discharge.

### 5.3. Advice for fire-fighters

As with all fires, wear positive pressure, self-contained breathing apparatus, (SCBA) with a full face piece and protective clothing. Persons without respiratory protection should leave area. Wear SCBA during clean-up immediately after fire. No smoking.

Fire fighters should wear appropriate protective equipment and self-contained breathing apparatus (SCBA) with a full face-piece operated in positive pressure mode.

Move containers away from fire area, if it can be done safely. Use water spray to cool fire-exposed containers. Do not allow firefighting water to enter drains or water courses. Fight fire with normal precautions from a reasonable distance.

ERG Guide No. 127

## Section 6. Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

Put on appropriate personal protective equipment (see section 8).

### 6.2. Environmental precautions

Do not allow spills to enter drains or waterways.

Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

### 6.3. Methods and material for containment and cleaning up

Prevent entry into sewers, water courses, basements or confined areas. Wash spillages into an effluent treatment plant or proceed as follows: Contain and collect spillage with non-combustible, absorbent material (e.g., sand, earth, vermiculite or diatomaceous earth) and place in container for disposal according to local regulations (see Section XIII). Use spark-proof tools and explosion-proof equipment. Dispose of via a licensed waste disposal contractor. Contaminated absorbent material may pose the same hazard as the spilled product. See Section I for emergency contact information and Section XIII for waste disposal.

## Section 7. Handling and storage

### 7.1. Precautions for safe handling

Handle containers carefully to prevent damage and spillage.

Put on appropriate personal protective equipment (see Section VIII). Eating, drinking and smoking should be prohibited in areas where this material is handled, stored and processed. Workers should wash hands and face before eating, drinking and smoking. Avoid contact with eyes, skin and clothing. Use only with adequate ventilation. Wear appropriate respirator, when ventilation is inadequate. Do not enter storage areas and confined spaces unless adequately ventilated. Use explosion-proof electrical equipment (ventilating, lighting and material handling). Use non-sparking tools. Take precautionary measures against electrostatic discharges. To avoid fire or explosion, dissipate static electricity during transfer by grounding and bonding containers and equipment before transferring material. Empty containers retain product residue and can be hazardous. Keep away from heat, sparks and flame.

Check section 2.2 (GHS Label Elements) for further details. - [Prevention]:



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## 7.2. Conditions for safe storage, including any incompatibilities

Store in accordance with local regulations. Store at room temperature. Store away from direct sunlight in a cool and well-ventilated place. Keep in properly labeled containers.

Incompatible materials: Reactive or incompatible with the following materials: oxidizing materials, acids and alkalis. Slightly reactive or incompatible with metals.

Check section 2.2 (GHS Label Elements) for further details. - [Storage]:

## 7.3. Specific end use(s)

No data available.

## Section 8. Exposure controls / personal protection

### 8.1. Control parameters

CAS No.	Ingredient	Exposure	
		Source	Value
0000064-17-5	Ethanol	OSHA	TWA 1000 ppm (1900 mg/m <sup>3</sup> )
		ACGIH	No Established Limit
		NIOSH	TWA 1000 ppm (1900 mg/m <sup>3</sup> )
0000077-86-1	Tris (hydroxymethyl) aminomethane	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit
0000593-84-0	Thiocyanic acid, compound with guanidine (1:1)	OSHA	No Established Limit
		ACGIH	No Established Limit
		NIOSH	No Established Limit

### 8.2. Exposure controls

<b>Respiratory</b>	If workers are exposed to concentrations above the exposure limit they must use the appropriate, certified respirators.
<b>Eyes</b>	Protective safety glasses recommended
<b>Skin</b>	Protective gloves recommended.
<b>Engineering Controls</b>	Provide adequate ventilation. Where reasonably practicable this should be achieved by the use of local exhaust ventilation and good general extraction. If these are not sufficient to maintain concentrations of particulates and any vapor below occupational exposure limits suitable respiratory protection must be worn.
<b>Other Work Practices</b>	Use good personal hygiene practices. Wash hands before eating, drinking, smoking or using toilet. Promptly remove soiled clothing and wash thoroughly before reuse.

Check section 2.2 (GHS Label Elements) for further details.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## Section 9. Physical and chemical properties

<b>Appearance</b>	Translucent blue low viscosity Liquid
<b>Odor</b>	Alcohol Odor
<b>Odor threshold</b>	Not determined
<b>pH</b>	7.9 – 8.3
<b>Melting point / freezing point</b>	Not Measured
<b>Initial boiling point and boiling range</b>	Not Measured
<b>Flash Point</b>	35 °C
<b>Evaporation rate (Ether = 1)</b>	Not Measured
<b>Flammability (solid, gas)</b>	Not Applicable
<b>Upper/lower flammability or explosive limits</b>	<b>Lower Explosive Limit:</b> Not Measured <b>Upper Explosive Limit:</b> Not Measured
<b>Vapor pressure (Pa)</b>	Not Measured
<b>Vapor Density</b>	Not Measured
<b>Relative Density</b>	1.06 – 1.10
<b>Solubility in Water</b>	Not Measured
<b>Partition coefficient n-octanol/water (Log Kow)</b>	Not Measured
<b>Auto-ignition temperature</b>	Not Measured
<b>Decomposition temperature</b>	Not Measured
<b>Viscosity (cSt)</b>	Not Measured
<b>9.2. Other information</b>	No other relevant information.

## Section 10. Stability and reactivity

### 10.1. Reactivity

Hazardous Polymerization will not occur.

### 10.2. Chemical stability

Stable under normal circumstances.

### 10.3. Possibility of hazardous reactions

No data available.

### 10.4. Conditions to avoid

Avoid all possible sources of ignition (spark or flame).

### 10.5. Incompatible materials

Reactive or incompatible with the following materials: oxidizing materials, acids and alkalis. Slightly reactive or incompatible with metals.

### 10.6. Hazardous decomposition products

Decomposition products may include carbon dioxide, carbon monoxide and nitrogen oxides.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## Section 11. Toxicological information

### Acute toxicity

Exposure to solvent vapor concentrations from the component solvents in excess of the stated occupational exposure limits may result in adverse health effects such as mucous membrane and respiratory system irritation and adverse effects on the kidneys, liver and central nervous system. Symptoms include headache, nausea, dizziness, fatigue, muscular weakness, drowsiness and in extreme cases, loss of consciousness.

Repeated or prolonged contact with the preparation may cause removal of natural fat from the skin resulting in dryness, irritation and possible non-allergic contact dermatitis. Solvents may also be absorbed through the skin. Splashes of liquid in the eyes may cause irritation and soreness with possible reversible damage.

Note: When no route specific LD50 data is available for an acute toxin, the converted acute toxicity point estimate was used in the calculation of the product's ATE (Acute Toxicity Estimate).

Ingredient	Oral LD50, mg/kg	Skin LD50, mg/kg	Inhalation Vapor LC50, mg/L/4hr	Inhalation Dust/Mist LC50, mg/L/4hr	Inhalation Gas LC50, ppm
Thiocyanic acid, compound with guanidine (1:1) - (593-84-0)	539.00, Rat - Category: 4	No data available	No data available	No data available	No data available
Ethanol - (64-17-5)	10,470.00, Rat - Category: NA	17,100.00, Rabbit - Category: NA	124.70, Rat - Category: NA	No data available	No data available
Tris (hydroxymethyl) aminomethane - (77-86-1)	No data available	No data available	No data available	No data available	No data available

### Carcinogen Data

CAS No.	Ingredient	Source	Value
0000064-17-5	Ethanol	OSHA	Regulated Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: Yes; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
		ACGIH	A3
0000077-86-1	Tris (hydroxymethyl) aminomethane	OSHA	Regulated Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
		ACGIH	No Established Limit
0000593-84-0	Thiocyanic acid, compound with guanidine (1:1)	OSHA	Regulated Carcinogen: No
		NTP	Known: No; Suspected: No
		IARC	Group 1: No; Group 2a: No; Group 2b: No; Group 3: No; Group 4: No;
		ACGIH	No Established Limit

Classification	Category	Hazard Description
Acute toxicity (oral)	4	Harmful if swallowed.
Acute toxicity (dermal)	5	May be harmful in contact with skin.
Acute toxicity (inhalation)	---	Not Applicable
Skin corrosion/irritation	3	Causes mild skin irritation.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

Serious eye damage/irritation	---	Not Applicable
Respiratory sensitization	---	Not Applicable
Skin sensitization	---	Not Applicable
Germ cell mutagenicity	---	Not Applicable
Carcinogenicity	---	Not Applicable
Reproductive toxicity	---	Not Applicable
STOT-single exposure	---	Not Applicable
STOT-repeated exposure	---	Not Applicable
Aspiration hazard	---	Not Applicable

## Section 12. Ecological information

### 12.1. Toxicity

No additional information provided for this product. See Section 3 for chemical specific data.

### Aquatic Ecotoxicity

Ingredient	96 hr LC50 fish, mg/l	48 hr EC50 crustacea, mg/l	ErC50 algae, mg/l
Thiocyanic acid, compound with guanidine (1:1) - (593-84-0)	Not Available	Not Available	Not Available
Ethanol - (64-17-5)	15,400.00, Lepomis macrochirus	> 10,000, Daphnia magna	17.921 (96 hr), Ulva pertusa
Tris (hydroxymethyl) aminomethane - (77-86-1)	Not Available	Not Available	Not Available

### 12.2. Persistence and degradability

There is no data available on the preparation itself.

### 12.3. Bioaccumulative potential

Not Measured

### 12.4. Mobility in soil

No data available.

### 12.5. Results of PBT and vPvB assessment

This product contains no PBT/vPvB chemicals.

### 12.6. Other adverse effects

No data available.

## Section 13. Disposal considerations

### 13.1. Waste treatment methods

Observe all federal, state and local regulations when disposing of this substance.





# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## Section 14. Transport information

	DOT (Domestic Surface Transportation)	IMO / IMDG (Ocean Transportation)	ICAO/IATA
<b>14.1. UN number</b>	UN1170	UN1170	UN1170
<b>14.2. UN proper shipping name</b>	UN1170, Ethanol solutions, 3, III	Ethanol solutions	Ethanol solutions
<b>14.3. Transport hazard class(es)</b>	<b>DOT Hazard Class: 3</b>	<b>IMDG: 3</b> <b>Sub Class: Not Applicable</b>	<b>Air Class: 3</b>
<b>14.4. Packing group</b>	III	III	III
<b>14.5. Environmental hazards</b>			
<b>IMDG</b>	Marine Pollutant: No;		
<b>14.6. Special precautions for user</b>	Not Applicable		

## Section 15. Regulatory information

<b>Regulatory Overview</b>	The regulatory data in Section 15 is not intended to be all-inclusive, only selected regulations are represented.
<b>Toxic Substance Control Act ( TSCA)</b>	All components of this material are either listed or exempt from listing on the TSCA Inventory.
<b>Note:</b>	Ethanol (Ethyl alcohol) has been placed on the California Proposition 65 Cancer and Developmental Lists when in alcoholic beverages.
<b>EPCRA 302 Extremely Hazardous:</b>	Hydrochloric acid
<b>EPCRA 313 Toxic Chemicals:</b>	Ethanol
<b>Proposition 65 - Carcinogens (&gt;0.0%):</b>	To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.
<b>Proposition 65 - Developmental Toxins (&gt;0.0%):</b>	To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.
<b>Proposition 65 - Female Repro Toxins (&gt;0.0%):</b>	To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.
<b>Proposition 65 - Male Repro Toxins (&gt;0.0%):</b>	To the best of our knowledge, there are no chemicals at levels which require reporting under this statute.
<b>Proposition 65 Label Warning:</b>	This product contains no chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.



# Safety Data Sheet

Product: SDNA-1000

SDS Revision Date: 06/15/2020

## Section 16. Other information

**SDS Revision Date** 06/15/2020

The information and recommendations contained herein are based upon data believed to be correct. However, no guarantee or warranty of any kind, expressed or implied, is made with respect to the information contained herein. We accept no responsibility and disclaim all liability for any harmful effects which may be caused by exposure to our products. Customers/users of this product must comply with all applicable health and safety laws, regulations, and orders.

The full text of the phrases appearing in section 3 is:

H225 Highly flammable liquid and vapor.

H302 Harmful if swallowed.

H312 Harmful in contact with skin.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H335 May cause respiratory irritation.

H412 Harmful to aquatic life with long lasting effects.

Disclaimer: The information presented herein is supplied as a guide to those who handle or use this product. Safe work practices must be employed when working with any materials. It is important that the end user makes a determination regarding the adequacy of the safety procedures employed during the use of this product.

End of Document