

KCNK5 Polyclonal Antibody

Catalog number: 14136-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

KCNK5; FLJ11035; K2p5.1;

TASK-2; TASK2



PC-3 cells were subjected to SDS PAGE followed by western blot with 14136-1-AP (KCNK5 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	PC-3 cells
Calculated KCNK5 MW:	499aa, 55kd
Observed KCNK5 MW:	72kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5355
GenBank accession number:	BC060793
Gene ID (NCBI):	8645
Full name:	Potassium channel, subfamily K, member 5

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.