

## MSL3 Polyclonal Antibody

Catalog number: 12977-1-AP

Size: 37 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

MSL3; DKFZp586J1822, Male

specific lethal 3 homolog,

MSL3, MSL3 like 1, MSL3L1



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 12977-1-AP(MSL3 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated MSL3 MW:	521aa,60 kDa
Observed MSL3 MW:	60 kDa
Positive WB detected in	Mouse skeletal muscle tissue, human stomach tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3592
GenBank accession number:	BC031210
Gene ID (NCBI):	10943
Full name:	Male-specific lethal 3 homolog (Drosophila)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.