

FCN3 Polyclonal Antibody

Catalog number: 11867-1-AP

Size: 25 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

FCN3; FCN3, FCNH, Ficolin 3,

HAKA1, Hakata antigen



human placenta tissue were subjected to SDS PAGE followed by western blot with 11867-1-AP(FCN3 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated FCN3 MW:	33 kDa
Observed FCN3 MW:	35kd+50kd
Positive WB detected in	Human placenta tissue, HepG2 cells, Y79 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2449
GenBank accession number:	BC020731
Gene ID (NCBI):	8547
Full name:	Ficolin (collagen/fibrinogen domain containing) 3 (Hakata antigen)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.