


GEB Universal tips for TWC thin wallcrown And for improved flexibility.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.5*3.5)mm	Packaging Spec.
PT0200-B-Y	200ul	Yellow				1000tips/bag,
PT0200-B-Y-LB	200ul	Yellow	Low Retention			20bags/case, 20000tips/case
Case dimension L*W*H: 595*325*310mm.			Gross Weight: 7.8kg.			
FT0200-B-N	200ul	Natural			Filtered	1000tips/bag,
FT0200-B-N-LB	200ul	Natural	Low Retention		Filtered	10bags/case,
FT0200-B-NS	200ul	Natural		Sterilized	Filtered	10000tips/case
FT0200-B-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 480*345*280mm.			Gross Weight: 4.5kg.			
PT0200-9-Y	200ul	Yellow				96tips/rack,
PT0200-9-Y-LB	200ul	Yellow	Low Retention			10racks/pack,
PT0200-9-YS	200ul	Yellow		Sterilized		5packs/case,
PT0200-9-YS-LB	200ul	Yellow	Low Retention	Sterilized		4800tips/case.
FT0200-R-NS	200ul	Natural		Sterilized	Filtered	50racks/case.
FT0200-R-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 500*265*450mm.			Gross Weight: 7.8kg.			

- ✧ 200μL with graduation marks at 10ul,50ul and 100ul
- ✧ 96 Tip/Hanging Tip Rack format / Bulk
- ✧ Polypropylene
- ✧ Presterile/Nonsterile
- ✧ Filtered (With aerosol-resistant filters) /NonFiltered
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Universal tips workstations.

Drawing: 200ul



Prepared by: Anne Oct.19, 2020 Approve: WeiMing Yu Oct.19,2020