

## EIF3K Polyclonal Antibody

Catalog number: 10640-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

EIF3K; ARG134, eIF 3 p25, eIF 3

p28, EIF3 p28, EIF3K, EIF3S12,

HSPC029, M9, MSTP001, PLAC

24, PLAC24, PRO1474, PTD001

### Applications

**Tested applications:** ELISA  
**Species specificity:** Human, Mouse, Rat; other species not tested.  
**Calculated EIF3K MW:** 28 kDa  
**Observed EIF3K MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Ag1000  
**GenBank accession number:** BC007559  
**Gene ID (NCBI):** 27335  
**Full name:** Eukaryotic translation initiation factor 3, subunit K

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.