



# Tube Revolvers

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

## Powerful Tube Revolver with Five Unique Interchangeable Rotisseries



### Product Features:

- Small footprint with easy operation
- Includes 5 unique interchangeable rotisseries
- Rotisseries can be mixed to accommodate a variety of tube sizes simultaneously
- Paddle angles are adjustable 0-90°
- Maintenance-free motor is powerful and quiet
- TR-02UA can be programmed to oscillate at every 40° in the rotation to simulate hand mixing with rotation and shaking

Both models include the following 5 rotisserie paddles



### Specifications:

Description	Tube Revolver with Fixed Speed	Tube Revolver with Variable Speed and Digital Display
Model	TR-01A	TR-02UA
Speed	25 rpm	10~40rpm
Speed Accuracy	N/A	±1rpm
Capacity	84x 0.5 / 0.8ml Tubes, 60x 1.5 / 2.0ml Tubes, 28x 5 / 7ml Tubes, 28x 10 / 15ml Tubes, 6x 50ml Tubes	
Display	N/A	LED
Included Accessories	TRA01 (2Pcs), TRA02 (2Pcs), TRA03 (2Pcs), TRA04 (2Pcs), TRA05 (1Pcs)	
Ambient Temperature	4°C~60°C	
Ambient Relative Humidity	≤80%	
Power Supply	AC110V, 60Hz	AC110V, 60Hz
Max. Power	4W	10W
Dimensions (DxWxH)	148x225x195mm (5.83x8.86x7.68")	148x260x195mm (5.83x10.24x7.68")
Net Weight	0.8kg (1.76lb)	1kg (2.2lb)
Protection Level	IP21	
Certifications	CE, RoHS, WEEE	CE, RoHS, WEEE, cCSAus

