

## UGCGL2 Polyclonal Antibody

Catalog number: 13420-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

UGCGL2; FLJ10873; FLJ11485;

HUGT2; MGC117360;

MGC150689; MGC87276



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13420-1-AP(UGCGL2 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Positive WB detected in:	HEK-293 cells
Calculated UGCGL2 MW:	1516aa,168kd
Observed UGCGL2 MW:	170kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag4189
GenBank accession number:	BC032302
Gene ID (NCBI):	55757
Full name:	UDP-glucose ceramide glucosyltransferase-like 2

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.