

STAMBP Polyclonal Antibody

Catalog number: 11346-1-AP


Size: 50 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

STAMBP; AMSH, STAM binding protein, STAMBP



Immunohistochemical of paraffin-embedded human prostate cancer using 11346-1-AP (STAMBP antibody) at dilution of 1:100 (under 10x lens)



mouse brain tissue were subjected to SDS PAGE followed by western blot with 11346-1-AP (STAMBP antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB, IHC
Cited applications:	WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Human
Calculated STAMBP MW:	48 kDa
Observed STAMBP MW:	48 kDa
Positive WB detected in	Mouse brain tissue, human placenta tissue
Positive IHC detected in	Human prostate cancer tissue
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1895
GenBank accession number:	BC007682
Gene ID (NCBI):	10617
Full name:	STAM binding protein

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.