

GRIPSTRONG® NITRILE

Surgical

Industrial

Exam

Non-Latex • Powder-Free • Textured



Sempermed®

We are people protection.®

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



GripStrong® Nitrile gloves provide superior gripping capabilities for multi-purpose protection at a great value, all without the risk of latex allergens. Advanced manufacturing technology gives GripStrong® Nitrile a softer fit and feel than traditional nitrile gloves allowing for maximum tactile sensitivity. GripStrong® is ideal for food, safety, janitorial/ sanitation, and laboratory markets.

FEATURES

- Manufactured with advanced technology
- Textured fingertips
- Supple touch
- Not made with natural rubber latex
- Powder-free
- ESD Tested and approved

BENEFITS

- Softer fit and feel than traditional nitrile gloves
- Improved dexterity and gripping capabilities
- Maximum tactile sensitivity and flexibility
- A latex alternative
- Eliminates glove powder aerosolization
- Safe for use with sensitive electronics

PRODUCT CHARACTERISTICS

- Synthetic Nitrile Rubber (NBR)
- Textured Fingertips
- Powder-Free
- Ambidextrous
- Beaded Cuff
- Single Use Only

Order #	Size	Gloves/Box*	Boxes/Case
GSNF102	S	100	10
GSNF103	M	100	10
GSNF104	L	100	10
GSNF105	XL	100	10

*By weight

All Sempermed gloves are single use only.



HOW WE GLOVE THE PLANET WITH GRIPSTRONG® NITRILE

- ✓ Inner boxes made with 80% recycled material
- ✓ Outer cartons made with 100% recycled material

Sempermed®

We are people protection.®

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558



Glove the Planet® is Sempermed USA Inc.'s commitment to reducing our environmental impact. By thinking ahead and consuming less we're able to do more to minimize our global footprint with every step we take. Our ethos of sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit www.GlovethePlanet.com.



TM/® Trademark/registered trademark of Sempermed USA, Inc. 06182016