



Laboratory Refrigerator/Freezer Combo For Sample Storage

Futura Silver Series combo refrigerator/freezer units for sample storage are ideal for laboratories with limited space. These stacked style combo units offer superior temperature performance and sample security to ensure safe storage at uniform temperatures.



Standard Product Features:

- Manual or Auto defrost freezer with cycle defrost refrigerator
- -13°F (-25°C) freezer operating temperature;
+39.2°F (4°C) refrigerator operating temperature
- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- White powder coat exterior (LHP-9-RFCSS-PH stainless steel exterior)
- All combination units are stacked and come with separate compressors, power cords and temperature controllers
- HFC-free refrigerant (R600a), foam insulation and packaging
- Keyed lock doors
- Magnetic door gaskets for positive seal
- Probe access ports (3/8") - one for refrigerator, one for freezer
- Pharmacy refrigerator/freezer toolkit
- 115V operation (requires two outlets)
- CDC approved power cords
- Remote alarm contacts
- Heavy duty wire shelves
- Leveling legs
- Warranties: 2 years parts and labor, 5 years compressor parts only.



**Microprocessor Temperature
Controllers**



Hydrocarbon refrigerant used in
all combo LHP models

LABRepCo
An Employee-Owned Company

Specifications and Ordering Information:



LHP-9-RFCSS

LHP-7-RFCG

Model		LHP-7-RFCG	LHP-7-RFCAG	LHP-9-RFCG	LHP-9-RFCSS	LHP-9-RFCSD
Weight lbs		135	135	200	200	200
Overall Height in.		53 1/8	53 1/8	66 3/4	66 3/4	66 3/4
Overall Width in.		23 5/8	23 3/4	23 3/4	23 3/4	23 3/4
Overall Depth in.		24 1/4	24 1/4	25 1/2	24 1/2	24 1/2
Interior Height in.	Refrig.	26 3/4	26 3/4	24 1/2	24 1/2	24 1/2
	Freezer	13 7/8	13 7/8	24 1/2	24	24 1/2
Interior Width in.	Refrig.	19 7/8	19 7/8	19 1/4	19 1/4	19 1/4
	Freezer	17 3/4	18	19	19	19
Interior Depth in.	Refrig.	16	16	17 1/4	17 1/2	17 1/4
	Freezer	16	9 1/8	17 1/4	17 1/4	17 1/4
Interior cu. ft.	Refrig.	5.2	5.2	4.6	4.5	4.6
	Freezer	1.7	1.3	4.2	4.2	4.2
# of Shelves:	Refrig.	3 adj	3 adj	1 fixed; 2 adj	1 fixed; 2 adj	1 fixed; 2 adj
	Freezer	2 adj	2 adj	2 fixed	2 fixed	2 fixed
Defrost Type	Refrig.	Cycle	Cycle	Cycle	Cycle	Cycle
	Freezer	Manual	Auto	Manual	Manual	Manual
Temperature:	Refrig.	+1°C to +10°C	+1°C to +10°C	+1°C to +10°C	+1°C to +10°C	+1°C to +10°C
	Freezer	-15°C to -25°C	-15°C to -25°C	-15°C to -25°C	-15°C to -25°C	-15°C to -25°C
Compressor Size HP	Refrig.	1/6	1/6	1/8	1/8	1/8
	Freezer	1/8	1/8	1/8	1/8	1/8
Plug (Receptacle)		5-15P (5-15R)	5-15P (5-15R)	5-15P (5-15R)	5-15P (5-15R)	5-15P (5-15R)
Total Amp Draw	Refrig.	1.3	1.3	1.74	1.74	1.74
	Freezer	1	1	1.5	1.5	5
Voltage		115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz



LHP-12-RFCAG

LHP-12-RFCSD

Model		LHP-12-RFCSD	LHP-12-RFCASD	LHP-12-RFCAG	LHP-12-RFCG
Weight lbs		219	219	246	246
Overall Height in.		80 1/4	80 1/4	80 1/4	80 1/4
Overall Width in.		23 5/8	23 3/4	23 3/4	23 5/8
Overall Depth in.		26 1/2	26 1/2	26 1/2	26 1/2
Interior Height in.	Refrig.	53 1/2	53 1/2	53 1/2	53 1/2
	Freezer	13 7/8	13 7/8	13 7/8	13 7/8
Interior Width in.	Refrig.	19 7/8	19 7/8	19 7/8	19 7/8
	Freezer	17 3/4	18	18	17 3/4
Interior Depth in.	Refrig.	18 3/8	18 3/8	18 3/8	18 3/8
	Freezer	16	9 1/8	9 1/8	16
Interior cu. ft.	Refrig.	10.5	10.5	10.5	10.5
	Freezer	1.7	1.3	1.3	1.7
# of Shelves:	Refrig.	6 adj	6 adj	6 adj	6 adj
	Freezer	2 adj	2 adj	2 adj	2 adj
Defrost Type	Refrig.	Cycle	Cycle	Cycle	Cycle
	Freezer	Manual	Auto	Auto	Manual
Temperature:	Refrig.	+1°C to +10°C	+1°C to +10°C	+1°C to +10°C	+1°C to +10°C
	Freezer	-15°C to -25°C	-15°C to -25°C	-15°C to -25°C	-15°C to -25°C
Compressor Size HP	Refrig.	1/3	1/3	1/3	1/3
	Freezer	1/8	1/8	1/8	1/8
Plug (Receptacle)		5-15P (5-15R)	5-15P (5-15R)	5-15P (5-15R)	5-15P (5-15R)
Total Amp Draw	Refrig.	2.3	2.3	2.3	2.3
	Freezer	1	1	1	1
Voltage		115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

Please Note: All dimensions are nominal. All specifications subject to change without notice.

Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15R

LABRepCo
An Employee-Owned Company