

BinLiner™ Stands



Micronova's **BinLiner™ Stands** are the perfect complement to our portfolio of BinLiner bags. There are two stand choices, rectangular or scissor style frame made of durable stainless steel. Both styles allow for simple, ergonomic side removal of liners.

The unique rectangular stand features polypropylene corner molds that connect the sides and function as the stand base. The base reduces the possibility of rusting directly on the floor surface and slides easily to help prevent scratching. Optional casters can be added for ease of mobility when relocating the bin stand. The scissor stand folds open and closes for storage when not in use.



(Left) Rectangular Stand
BLSSH-1624 with 40x48 in liner



(Right)
Scissor stand
BLSSH-40 with
40x48 in liner

SPECIFICATIONS

- Two BinLiner Stand options: Rectangular stainless steel frame with polypropylene corners or scissor style stainless steel frame
- Securely holds a range of bin liners for material disposal, garment collection and raw material or product storage
- Ergonomic side removal simplified when removing liners
- Corner legs are removable for easy disassembly and cleaning. Rectangular stand available with optional casters, as special order
- Available stands hold bags from 4 to 200 L (1 to 50-gallon) capacities.
- Resistant to common cleaning solutions and disinfectants
- Autoclavable
- Good chemical resistance to disinfectants

Rectangular Stand
BLSSH-1624



Manufactured for use in Controlled Environments

- ISO 5 thru 9 Cleanrooms
- ISO 9001 Certified Manufacturing Facility

Materials of Construction:

Rectangular Stand: 304 stainless steel ¼" rods -electropolished.
Polypropylene corner connectors
Scissor Stand: Stainless Steel

Dimensions:

BLSSH-46

Top & Bottom: 5 x 5 in (12.7 x 12.7 cm)
Height: 11 in (28 cm)

Use with:
10 x 16 in size bags

BLSSH-168

Top & Bottom: 16 x 8 in (40.6 x 20.3 cm)
Height: 28 in (71 cm)

Use with:
24 x 36 in and
25 x 36 in size bags

BLSSH-1515

Top & Bottom: 15 x 15 in (38 x 38 cm)
Height: 28 in (71 cm)

Use with:
34 x 40 in size bags

BLSSH-1624*


Top & Bottom: 16 x 25 in (40.6 x 63.5 cm)
Height: 34 in (86.6 cm)

Use with:
40 x 48 in and
40 x 60 in size bags*

Scissor Stand:
BLSSH-40

Use with:
40 x 48 in size bags

* For 6mil 40 x 48 in and 40 x 60 in bags -
use frame with white tape as top edge



Micronova Manufacturing
3431 W. Lomita Boulevard
Torrance, CA 90505-5010 USA

Tel: 310/784-6990
Fax: 310/784-6980

SalesInfo@micronova-mfg.com
www.micronova-mfg.com