

SAVE 10% on our New Wiley® Mini Mills!



Single Speed



Variable Speed

- Power: 115V 60Hz or 230V 50Hz
- Speed: 1200 rpm free spindle speed
- Chamber Size: 1.575" Diameter
- Stationary Blades: 2
- Rotating Blades: 4
- Max Sample Input Size: .5" Diameter
- Continuous Torque: 1.4 Nm
- Exterior Material: Powder Coated Steel

- Power: 115V 60Hz or 230V 50Hz
- Speed: 1-2400 rpm free spindle speed
- Chamber Size: 1.575" Diameter
- Stationary Blades: 2
- Rotating Blades: 4
- Max Sample Input Size: .5" Diameter
- Continuous Torque: 1.4 Nm
- Exterior Material: Powder Coated Steel

The Mini Mills are all steel Benchtop Cutting Mills featuring a 1.575-inch diameter sealed cutting chamber with two adjustable stationary blades and a rotating cutting head with four cutting edges that revolve at high speed to produce a shearing action that minimizes heat build-up.

The hard, steel cutting edges on the knives allow milling of a wider range of materials, including plastics, animal feed, bones, cables, cardboard, electronic scrap, feed pellets, foils, food, leather, material mixtures, non-ferrous metals, paper, pharmaceutical products, refuse-derived fuels, rubber, waste, and wood.

Thomas No.	Description	List Price	Sale Price
23A00D287	Wiley® Mini Mill, Single-Speed, 115V 60HZ (replaces 1171H10)	\$6,650.00	\$5,985.00
23A00D288	Wiley® Mini Mill, Single-Speed, 230V 50HZ	\$7,736.63	\$6,962.97
23A00D289	Wiley® Mini Mill, Variable-Speed, 115V 60HZ (replaces 1173U48)	\$10,341.00	\$9,306.90
23A00D290	Wiley® Mini Mill, Variable-Speed, 230V 50HZ (replaces 1173U48)	\$10,942.70	\$9,848.43

**Use
PROMO Code:
0801
at checkout**

Special Pricing Valid September 1-30, 2023!