

## TSGA10 Polyclonal Antibody

Catalog number: 12593-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TSGA10; CEP4L



mouse testis tissue were subjected to SDS PAGE followed by western blot with 12593-1-AP(TSGA10 Antibody) at dilution of 1:1500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse testis tissue, mouse brain tissue
Calculated TSGA10 MW:	698aa,81kd
Observed TSGA10 MW:	81kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3311
GenBank accession number:	BC028366
Gene ID (NCBI):	80705
Full name:	Testis specific, 10

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.