# **Compact Glove Boxes**

Compact design with transfer chamber and flat side access door



## SG828/835/848/860

## A completely sealed glove box for PRODUCT AND PERSONNEL PROTECTION





**SG-828** 28" W x 23" D x 29" H 710 x 580 x 740 mm

**SG-860** 60" W x 29" D x 32" H 1524 x 740 x 810 mm

All models feature high visibility interior, cleanliness, safety, complete containment, portability and a draft-free atmopshere

All models are completed systems, nothing needs to be added except gas of choice: nitrogen, argon or other inert gases Installed purging valves are ideal for lowering oxygen and himidity levels

Very useful when working with toxic substances, asbestos, fibers, sewage residue and harmful liquid vapors

### Features

- Constructed with an optically clear material, developed for isolation and/or containment
- Transfer chamber with purging valves and vacuum gauge
- Flat side access door for easy introduction of larger equipment
- Four (4) key cock valves for purging: 2 on the main body, 2 on the transfer chamber
- Hospital grade multiple electrical outlet strip
- Pressure relief valve with small HEPA filer
- White ambidextrous hypalon gloves

#### Specifications

Model	SG-828	SG-835	SG-848	SG-860
<b>RECOMMENDED OPERATIONAL PRE</b>	SSURE			
For containment purposes	-0.5 of water column (0.93 torr)			
For isolation purposes	0.5" of water column (0.93 torr)			
Main chamber				
Maximum pressure	+6" of water column (11.2 torr)			
Maximum vacuum	-6" of water column (11.2 torr)			
Transfer chamber				
Maximum pressure	Not engineered to support positive pressure			
Maximum vacuum	-26" of Hg (660 torr)			
Inside dimensions (W x D x H)	28" x 23" x 29" 710 x 580 x 740 mm	35" x 29"x 30" 890 x 740 x 760 mm	48" x 29" x 32" 1220 x 740 x 810 mm	60" x 29" x 32" 1524 x 740 x 810 mm
Outside dimensions (Includes transfer chamber 12" long)	43" x 24" x 31" 1100 x 610 x 790 mm	49" x 30" x 31" 1250 x 760 x 790 mm	63" x 31" x 35" 1600 x 790 x 860 mm	76" x 32" x 35" 1930 x 812 x 890 mm
Access door opening	15.5" x 22" h 400 x 560 mm	21.5" x 22" h 546 x 560 mm	21.5" x 24" h 546 x 610 mm	21.5" x 24" h 546 x 610 mm