



## Safety Data Sheet

### SECTION 1: Identification

#### 1.1. Product Identifier

**Trade Name or Designation:** VeriSpec™ Mixed Anion Standard 7

2 ppm F<sup>-</sup>, 10 ppm PO<sub>4</sub><sup>3-</sup>, NO<sub>3</sub><sup>-</sup>, Br<sup>-</sup>, SO<sub>4</sub><sup>2-</sup>, 5ppm Cl<sup>-</sup>, NO<sub>2</sub><sup>-</sup>, Manufactured and Tested in an IS

**Product Number:** RV010738

**Other Identifying Product Numbers:** RV010738-100N

#### 1.2. Recommended Use and Restrictions on Use

Calibration Standard

#### 1.3. Details of the Supplier of the Safety Data Sheet

**Company:** Ricca Chemical Company

**Address:** 448 West Fork Drive  
Arlington, TX 76012 USA

**Telephone:** 888-467-4222

#### 1.4. Emergency Telephone Number (24 hr)

CHEMTREC (USA) 800-424-9300  
CHEMTREC (International) 1+ 703-527-3887

### SECTION 2: Hazard(s) Identification

#### 2.1. Classification of the Substance or Mixture (in accordance with OSHA HCS 29 CFR 1910.1200)

*For the full text of the Hazard and Precautionary Statements listed below, see Section 16.*

This product is not categorized as hazardous in any GHS hazard class.

#### 2.2. GHS Label Elements

**Pictograms:** None required.

**Signal Word:** None required.

**Hazard Statements:** None required.

**Precautionary Statements:** None required.

## Safety Data Sheet

### 2.3. WHMIS Classification

WHMIS classification is not included based on the recommended option (Option 4) found in the Canada Gazette Part II, Vol. 149, No.3, page 458

### 2.4. Hazards not Otherwise Classified or Covered by GHS

Data not available.

## SECTION 3: Composition / Information on Ingredients

### 3.1. Components of Substance or Mixture

Chemical Name	Formula	Molecular Weight	CAS Number	Weight%
Water	H <sub>2</sub> O	18.01 g/mol	7732-18-5	1.00%
Sodium Fluoride	NaF	41.98 g/mol	7681-49-4	0.00%
Sodium Nitrite	NaNO <sub>2</sub>	68.99 g/mol	7632-00-0	0.00%
Sodium Chloride	NaCl	58.44 g/mol	7647-14-5	0.00%
Sodium Nitrate	NaNO <sub>3</sub>	84.99 g/mol	7631-99-4	0.00%
Ammonium Sulfate	(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	132.13 g/mol	7783-20-2	0.00%
Potassium Bromide	KBr	119.00 g/mol	7758-02-3	0.00%
Ammonium Dihydrogen Phosphate	NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub>	97.99 g/mol	7722-76-1	0.01%

## SECTION 4: First-Aid Measures

### 4.1. General First Aid Information

**Eye Contact:** No action required to be taken. If necessary, rinse eyes with water.

**Inhalation:** Not expected to require first aid. If necessary, remove to fresh air.

**Skin Contact:** No action required to be taken. If necessary, wash areas of contact with water.

**Ingestion:** No action required to be taken. If necessary, dilute with water.

### 4.2. Most Important Symptoms and Effects, Acute and Delayed

May cause mild irritation to areas of contact.

### 4.3. Medical Attention or Special Treatment Needed

Not expected to require special treatment.

## SECTION 5: Fire-Fighting Measures

### 5.1. Extinguishing Media

Not considered to be a fire or explosion hazard.



## Safety Data Sheet

### 5.2. Specific Hazards Arising from the Substance or Mixture

Not considered to be a fire or explosion hazard.

### 5.3. Special Protective Equipment for Firefighters

Wear protective clothing and NIOSH-approved breathing equipment appropriate for the surrounding fire.

## SECTION 6: Accidental Release Measures

### 6.1. Personal Precautions, Protective Equipment and Emergency Procedures

Wear appropriate PPE for the size and nature of the spill. As a general rule, wear safety glasses and gloves.

### 6.2. Cleanup and Containment Methods and Materials

Absorb with suitable material and dispose of in accordance with local regulations.

## SECTION 7: Handling and Storage

### 7.1. Precautions for Safe Handling and Storage Conditions

Protect from freezing and physical damage.

## SECTION 8: Exposure Controls / Personal Protection

### 8.2. Exposure Controls

**Engineering Controls:** No specific controls are needed. Normal room ventilation is adequate.

**Respiratory Protection:** No specific controls are needed. Normal room ventilation is adequate.

**Skin Protection:** No specific controls are needed.

**Eye Protection:** No specific controls are needed.

### 8.3. Personal Protective Equipment

As a general rule, wear safety glasses and gloves.



## Safety Data Sheet

### SECTION 9: Physical and Chemical Properties

#### 9.1. Basic Physical and Chemical Properties

**Appearance:** Data not available.

**Physical State:** Liquid

**Odor:** Data not available.

**Odor Threshold:** Data not available.

**pH:** Data not available.

**Melting/Freezing Point:** Data not available.

**Initial Boiling Point /Range:** Data not available.

**Flash Point:** Data not available.

**Evaporation Rate:** Data not available.

**Flammability:** Data not available.

**Flammability/Explosive Limits:** Data not available.

**Vapor Pressure:** Data not available.

**Vapor Density:** Data not available.

**Relative Density:** 1.07

**Solubility:** Data not available.

**Partition Coefficient (n-Octanol/Water):** Data not available.

**Auto-Ignition Temperature:** Data not available.

**Decomposition Temperature:** Data not available.

**Viscosity:** Data not available.

**Explosive Properties:** Data not available.

**Oxidizing Properties:** Data not available.

### SECTION 10: Stability and Reactivity

#### 10.1. Reactivity and Chemical Stability

Stable under normal conditions of use and storage.

#### 10.2. Possibility of Hazardous Reactions

Data not available.

#### 10.3. Conditions to Avoid and Incompatible Materials

Protect from freezing and physical damage.

#### 10.4. Hazardous Decomposition Products

May emit irritating fumes when heated to decomposition.



## Safety Data Sheet

### SECTION 11: Toxicological Information

#### 11.1. Information on Toxicological Effects

**Acute Toxicity - Oral Exposure:**

Not applicable.

**Acute Toxicity - Dermal Exposure:**

Not applicable.

**Acute Toxicity - Inhalation Exposure:**

Not applicable.

**Acute Toxicity - Other Information:**

Data not available.

**Skin Corrosion and Irritation:**

Not applicable.

**Serious Eye Damage and Irritation:**

Not applicable.

**Respiratory Sensitization:**

Not applicable.

**Skin Sensitization:**

Not applicable.

**Germ Cell Mutagenicity:**

Not applicable.

**Carcinogenicity:**

Not applicable.

**Reproductive Toxicity:**

Not applicable.

**Specific Target Organ Toxicity from Single Exposure:**

Not applicable.

**Specific Target Organ Toxicity from Repeated Exposure:**

Not applicable.

**Aspiration Hazard:**

Not applicable.

**Additional Toxicology Information:**

Data not available.

## **Safety Data Sheet**

### **SECTION 12: Ecological Information**

#### **12.1. Ecotoxicity**

Not applicable.

#### **12.2. Persistence and Degradability**

Data not available.

#### **12.3. Bioaccumulative Potential**

Data not available.

#### **12.4. Mobility in Soil**

Data not available.

#### **12.5. Other Adverse Ecological Effects**

Data not available.

### **SECTION 13: Disposal Considerations**

#### **13.1. Waste Treatment Methods**

Data not available.

### **SECTION 14: Transportation Information**

#### **14.1. Transportation by Land - Department of Transportation (DOT, United States of America)**

Not regulated according to DOT Regulations.



## Safety Data Sheet

Not regulated according to IATA Regulations.

### SECTION 15: Regulatory Information

#### 15.1. Occupational Safety and Health Administration (OSHA) Hazards

Not listed.

#### 15.2. Superfund Amendments and Reauthorization Act (SARA) 302 Extremely Hazardous Substances

Not listed.

#### 15.3. Superfund Amendments and Reauthorization Act (SARA) 311/312 Hazardous Chemicals

Sodium Nitrite (CAS # 7632-00-0): 100 lb final RQ; 45.4 kg final RQ

Sodium Fluoride (CAS # 7681-49-4): 1000 lb final RQ; 454 kg final RQ

#### 15.4. Superfund Amendments and Reauthorization Act (SARA) 313 Toxic Release Inventory (TRI)

Sodium Nitrate (CAS # 7631-99-4): 1.0 % de minimis concentration (reportable only when in aqueous solution, Chemical Category N511)

Sodium Nitrite (CAS # 7632-00-0): 1.0 % de minimis concentration

Ammonium Dihydrogen Phosphate (CAS # 7722-76-1): 1.0 % de minimis concentration (10% of total aqueous Ammonia is reportable under this listing)

Ammonium Sulfate (CAS # 7783-20-2): 1.0 % de minimis concentration (10% of total aqueous Ammonia is reportable under this listing)

#### 15.5. Massachusetts Right-to-Know Substance List

Sodium Nitrate (CAS # 7631-99-4): Present

Sodium Nitrite (CAS # 7632-00-0): Present

Sodium Fluoride (CAS # 7681-49-4): Present

Ammonium Sulfate (CAS # 7783-20-2): Present

#### 15.6. Pennsylvania Right-to-Know Hazardous Substances

Sodium Nitrate (CAS # 7631-99-4): Present

Sodium Nitrite (CAS # 7632-00-0): Environmental hazard

Sodium Nitrite (CAS # 7632-00-0): Present

Sodium Fluoride (CAS # 7681-49-4): Environmental hazard

Sodium Fluoride (CAS # 7681-49-4): Present

Water (CAS # 7732-18-5): Present

Ammonium Sulfate (CAS # 7783-20-2): Environmental hazard

Ammonium Sulfate (CAS # 7783-20-2): Present

## Safety Data Sheet

### 15.7. New Jersey Worker and Community Right-to-Know Components

Sodium Nitrate (CAS # 7631-99-4): sn 3722  
Sodium Nitrate (CAS # 7631-99-4): SN 3722 500 lb TPQ (water dissociable, Category Code N511)  
Sodium Nitrite (CAS # 7632-00-0): sn 2258  
Sodium Nitrite (CAS # 7632-00-0): SN 2258 500 lb TPQ  
Sodium Fluoride (CAS # 7681-49-4): sn 0936  
Sodium Fluoride (CAS # 7681-49-4): sn 1699  
Sodium Fluoride (CAS # 7681-49-4): teratogen

### 15.8. California Proposition 65

Not listed.

### 15.9. Canada Domestic Substances List / Non-Domestic Substances List (DSL/NDSL)

Sodium Nitrate (CAS # 7631-99-4): Present (DSL)  
Sodium Nitrite (CAS # 7632-00-0): Present (DSL)  
Sodium Chloride (CAS # 7647-14-5): Present (DSL)  
Sodium Fluoride (CAS # 7681-49-4): Present (DSL)  
Ammonium Dihydrogen Phosphate (CAS # 7722-76-1): Present (DSL)  
Water (CAS # 7732-18-5): Present (DSL)  
Potassium Bromide (CAS # 7758-02-3): Present (DSL)  
Ammonium Sulfate (CAS # 7783-20-2): Present (DSL)

### 15.10. United States of America Toxic Substances Control Act (TSCA) List

Sodium Nitrate (CAS # 7631-99-4): Present  
Sodium Nitrite (CAS # 7632-00-0): Present  
Sodium Chloride (CAS # 7647-14-5): Present  
Sodium Fluoride (CAS # 7681-49-4): Present  
Ammonium Dihydrogen Phosphate (CAS # 7722-76-1): Present  
Water (CAS # 7732-18-5): Present  
Potassium Bromide (CAS # 7758-02-3): Present  
Ammonium Sulfate (CAS # 7783-20-2): Present

### 15.11. European Inventory of Existing Commercial Chemical Substances (EINECS), European List of Notified Chemical Substances (ELINCS), and No Longer Polymers (NLP)

Not listed.

## SECTION 16: Other Information

### 16.1. Full Text of Hazard Statements and Precautionary Statements



## Safety Data Sheet

### 16.2. Miscellaneous Hazard Classes

**Canadian Carcinogenicity Hazard Class:** Not Applicable.


**Physical Hazards Not Otherwise Classified (PHNOC):** Not Applicable.

**Health Hazards Not Otherwise Classified (HHNOC):** Not Applicable.

**Biohazardous Infectious Materials Hazard Class:** Not Applicable.

### 16.3. National Fire Protection Association (NFPA) Rating

**Health:** 0  
**Flammability:** 0  
**Reactivity:** 0  
**Special Hazard:**



## Safety Data Sheet

### 16.4. Document Revision

**Last Revision Date:** 9/15/2016

### DISCLAIMER

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and RICCA CHEMICAL COMPANY assumes no legal responsibility or liability whatsoever resulting from its use.