


Clarity Specialty Drug Testing Dip Cards





Easy-to-use dip cards



Custom multi drug cups and panels available with specialty drugs



Forensic use only

Drug Tested	Abbreviation	Cut-off
6-Monoacetylmorphine	6-MAM	10 ng/mL
Cathinones (Bath salts)*	MDPV	1000 ng/mL
Ethyl Glucuronide*	ETG	500 ng/mL
Fentanyl*	FTN	20 ng/mL
K2/Spice*	K2	50 ng/mL
Ketamine*	KET	1000 ng/mL
Lyserguc acid diethylamide*	LSD	20 ng/mL
Methaqualone*	MQL	300 ng/mL
Methylphenidate*	MPD	1000 ng/mL (300 ng/mL available)
Tenamfetamine*	MDA	500 ng/mL
Tramadol*	TRA	200 ng/mL
Urine Alcohol*	ALC	0.02%/0.04%/0.08%/0.30%
Zolpidem*	ZOL	50 ng/mL


^{*}Forensic Use Only

For Forensic IVD use, Clinical Trials, Pre-employement, and Employee screening only.

Read the drug test results at 5 minutes, do not interpret the results after 10 minutes.







Clear Vision - Clear Results

Clarity Brand Rapid Diagnostics