

Synaptotagmin-3 Polyclonal Antibody

Catalog number: 12726-1-AP

Size: 26 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

SYT3; DKFZp761O132,

Synaptotagmin 3,

synaptotagmin III,

Synaptotagmin-3, SYT3, SytIII



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 12726-1-AP(SYT3 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated Synaptotagmin-3	590aa, 63kd
MW:	
Observed Synaptotagmin-3	63kd
MW:	
Positive WB detected in	Mouse thymus tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3425
GenBank accession number:	BC028379
Gene ID (NCBI):	84258
Full name:	Synaptotagmin III

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.