

## CHRNA5 Polyclonal Antibody

Catalog number: 13516-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CHRNA5; CHRNA5, LNCR2,  
NACHRA5



Jurkat cells were subjected to SDS PAGE followed by western blot with 13516-1-AP(CHRNA5 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated CHRNA5 MW:	468aa,53 kDa
Observed CHRNA5 MW:	53 kDa
Positive WB detected in	Jurkat cells, A431 cells, A549 cells, mouse bladder tissue, mouse stomach tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4447
GenBank accession number:	BC033639
Gene ID (NCBI):	1138
Full name:	Cholinergic receptor, nicotinic, alpha 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.