

APBB2 Polyclonal Antibody

Catalog number: 13177-1-AP

Size: 47 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

APBB2; APBB2,

DKFZp434E033, FE65L, FE65L1,

Protein Fe65 like 1



human brain tissue were subjected to SDS PAGE followed by western blot with 13177-1-AP(APBB2 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated APBB2 MW:	83kd
Observed APBB2 MW:	83-95 kDa
Positive WB detected in	Human brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3844
GenBank accession number:	BC027946
Gene ID (NCBI):	323
Full name:	Amyloid beta (A4) precursor protein-binding, family B, member 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.