

RTKN Polyclonal Antibody

Catalog number: 12762-1-AP


Size: 43 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

RTKN; rhotekin, RTKN, RTKN1



Immunohistochemical of paraffin-embedded human prostate cancer using 12762-1-AP (RTKN antibody) at dilution of 1:100 (under 10x lens)



A431 cells were subjected to SDS PAGE followed by western blot with 12762-1-AP (RTKN antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated RTKN MW:	550aa, 61 kDa
Observed RTKN MW:	63kd
Positive WB detected in	A431 cells, mouse brain tissue, rat brain tissue
Positive IHC detected in	Human prostate cancer tissue
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3265
GenBank accession number:	BC017727
Gene ID (NCBI):	6242
Full name:	Rhotekin

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.