

The Reservoir You Need to Fill:

New and Sustainable Options



Reduce & Reuse

- Product core is made from reprocessed material, accounting for 50% of its composition
- Clearly marked with PS6 for recycling
- Domestically manufactured, reducing carbon footprint

Cooler Accessory

- Reusable Reagent Reservoir Cooler is a custom-fit companion for VistaLab 100mL Reagent Reservoirs
- Manufactured with non-toxic gel encased in heavy duty PVC

Small Footprint Reservoir



Small Footprint 50mL Reservoirs use 9% less plastic than standard 50mL Reservoirs




Mindful Packaging & Logistics

- Weekly pallets scheduled to minimize freight impact, reducing carbon emissions




Organizations We Support



To support sustainability efforts, the inner core is comprised of reprocessed resin

Uniquely treated to provide maximum reagent recovery, up to 67% better than comparable reservoirs!

The lowest residual volume in the industry



- 1 Trough-within-a-trough design to recover the "last few drops"
- 2 Pour-off spouts at every corner
- 3 Well-marked molded graduations on the inside wall (all sizes except 100mL)

Part No.	Product Description	Qty/Case
1188U07	10mL Reagent Reservoir, Bulk, Sterile	200
1226Y96	25mL Reagent Reservoir, Bulk, Sterile	200
1226Y98	50mL Reagent Reservoir, Bulk, Sterile	200
1226Y94	100mL Reagent Reservoir, Bulk, Sterile	200
1226Z02	100mL Reagent Reservoir Cooler	1
26A00F097	NEW 50mL Small Footprint Reagent Reservoir, Bulk, Sterile	100
26A00F100	NEW 50mL Small Footprint Reagent Reservoir, Bulk, Non-Sterile	100

*Other packaging options also available



CELLTREAT Scientific Products, LLC is a Committed Sustainable Supplier



833.544.SHIP (7447)
 CustomerService@ThomasSci.com
 www.ThomasSci.com

