

## TRIM37 Polyclonal Antibody

Catalog number: 13037-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TRIM37; TRIM37, POB1, TEF3,

MUL, tripartite motif

containing 37



PC-3 cells were subjected to SDS PAGE followed by western blot with 13037-1-AP (TRIM37 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated TRIM37 MW:	964aa, 108 kDa
Observed TRIM37 MW:	89kd
Positive WB detected in	PC-3 cells, mouse brain tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3690
GenBank accession number:	BC036012
Gene ID (NCBI):	4591
Full name:	Tripartite motif-containing 37

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.