

Perfect for cleanroom environments and strict sanitation protocols

These pre-cut strips are sterilized using a Vapor Hydrogen Peroxide process, adding an extra layer of protection to Parafilm's waterproof and moisture resistant properties

Ideal for protecting samples and preventing contamination or evaporation in sensitive environments

- Available in 2 pre-cut sizes
- Made, packaged, and sterilized in USA
- Non-toxic, odorless, and semi-transparent
- Creates a tight moisture barrier that holds a strong seal and doesn't release from glassware
- Prevents moisture and volume loss due to material composition and tight seal
- The flexible and self-sealing material allows it to cling around irregular shapes and surfaces
- Resistant up to 48 hours against saline solutions, inorganic acids, and alkaline solutions

STERILE | VH202


Case of 50 Packs

Item No.	Description	L x W x H in	L x W x H mm
234526ASTL	Parafilm® M, Sterile, 2" x 8" Strip, 10/PK	2 x 8	51 x 203
234526BSTL	Parafilm® M, Sterile, 4" x 6" Strip, 5/PK	4 x 6	101 x 152



Color sterilization indicator – turns pink when sterilized

Not sterile



Sterile



Sterile

