

OVOL1 Polyclonal Antibody

Catalog number: 14082-1-AP

Size: 38 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

OVOL1; HOVO1, ovo like

1(Drosophila), OVOL1



human brain tissue were subjected to SDS PAGE followed by western blot with 14082-1-AP(OVOL1 antibody) at dilution of 1:400



IP Result of anti-OVOL1 (IP:14082-1-AP, 4ug; Detection:14082-1-AP 1:500) with mouse brain tissue lysate 4000ug.

Applications

Tested applications:	ELISA, WB, IP
Cited applications:	WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Human
Calculated OVOL1 MW:	30 kDa
Observed OVOL1 MW:	30kd
Positive WB detected in	Human brain tissue, mouse pancreas tissue
Positive IP detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000 IP: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag5210
GenBank accession number:	BC059408
Gene ID (NCBI):	5017
Full name:	Ovo-like 1(Drosophila)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.