

LAB SAVINGS

Locking Microcentrifuge Tubes

0.5mL, 1.5mL, 2.0mL & 5mL Midi Tube

Globe Scientific locking microcentrifuge tubes are engineered with a secure locking system that provides the security of a screw cap microtube with the simplicity of a snap cap tube.

Secure – Maintains sample integrity with locking cap which prevents accidental cap opening and evaporation during centrifugation

Versatile – Single use tubes for laboratory procedures including centrifugation, boiling, incubation, freezing and sample storage

Quality – Tubes are made from optically clear USP Class VI medical grade polypropylene

Certified – Manufactured in a class 100,000 cleanroom and are certified to be non-pyrogenic, free from RNase, DNase, human DNA, PCR Inhibitors, ATP and heavy-metals

Ease of Use – Tubes feature graduations, writing area, pierceable cap, are autoclavable, and are packaged in zip, re-closable bags

**Buy 2
Get 1
FREE**



Thomas #	Item #	Description	Unit	Centrifuge Max. RCF
1169J23	111552LK	0.5mL, polypropylene (PP), natural,	500	30,000 x g
1169J24	111562LK	1.5mL, polypropylene (PP), natural	500	30,000 x g
1169J25	111572LK	2.0mL, polypropylene (PP), natural	500	30,000 x g
1169J26	111582LK	5.0mL, polypropylene (PP), natural	100	21,000 x g

BUY 2 CASES of Locking Microcentrifuge Tubes between April 1 and June 30, 2025, and GET 1 CASE FREE!

Limit 3 free cases. Mix & match. Free item of equal or lesser value. Offer good while supplies last. Valid for shipment within the contiguous US only.

For FREE Product(s), Please Complete Redemption Information on Reverse Side!

SPRING LAB SAVINGS



**Buy 2
Get 1
FREE**

ACCOUNT INFORMATION

First Name

Last Name

Date

Organization

Account No.

Shipping Address

Rm./Bldg./Dept.

City

State

Zip Code

Phone

Email



**MY FREE Product
Selection(s)**

Allow 3-4 weeks for delivery of
free products.

Valid 4/01/25 - 6/30/25

Complete the account information and
submit, along with Thomas Scientific proof of
purchase, to:

Promotions@thomassci.com

Reference **LOCKING2/1** when submitting redemption.

*Disclaimer: Limit 3 free cases. Mix & match. Free item of equal or
lesser value. Offer good while supplies last. Valid for shipment
within the contiguous US only. Redemptions must be received
no later than July 15, 2025. Pricing and offers subject to change
without notice.*

