



Graphing

Note: labeling is forbidden on the Freezing Bags body;
a small label can be placed on the surface of the inner PE bag
and a product label on the outer PE bag

Item No.	Qty	Description	Inventory Code	Material
D1	1	50mL Freezing Bags(0.13)	ACAX0022	EVA
A1	2	0.203"ID/0.280"OD(5.15mm/7.1mm)EVA	EABX0083	EVA
B1	1	2.95*4.15mm EVA Tubing,Length 18cm	DDAX0063	EVA
C1	1	Freeze tube tubing assembly	DDAX0060	PVC/Silicone/PE
F1	1	Light-curing glue	ADHX0001	PAA
F2	2	PE Bag	ADCX0011	PE
F3	1	2# Carton,,5pcs/pack, 20 packs/box	ADEX0001	Corrugated Paper

Customer Approval		LePure Biotech Co., Ltd.			
Company:		Drawing No. CR0050-1E			
Signature:		Title: SU Freezing Bags			
Date:		Designed By:	Date:	Checked By:	Date:
		Jiajia Guo	2023.03.17	ShuaikeLi	2023.03.17
Design Approval		Edition Rev03			
the manufacturing of the system(s) based on this design. By signing, you also assume full		We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property			