

## Dysbindin Polyclonal Antibody

Catalog number: 11132-1-AP

Size: 42 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DTNBP1; DBND, DTNBP1,


DTNBP1a, DTNBP1b,

Dysbindin, Dysbindin 1,

dystrobrevin binding protein 1,

HPS7, HPS7 protein, My031,

SDY



mouse brain tissue were subjected to SDS PAGE followed by western blot with 11132-1-AP(DTNBP1 antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated Dysbindin MW:	351aa,39 kDa
Observed Dysbindin MW:	48kd,36kd
Positive WB detected in	Mouse brain tissue, human testis tissue, mouse testis tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag1613
GenBank accession number:	BC011912
Gene ID (NCBI):	84062
Full name:	Dystrobrevin binding protein 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.