

Reagent Reservoirs

Choose from four styles:
25 mL Divided, 25 mL, 50 mL and 100 mL

- "V" bottom design allows for maximum recovery of reagents
- Graduated with convenient pouring spouts at all four corners
- Extra wide base for increased stability
- Unique "Nesting" feature makes separation of bulk packed reservoirs easy



25 mL divided reservoirs are divided into 3 separate 5 mL channels; great for low volume pipetting. Sterile and certified version are certified RNase and DNase-free and non-pyrogenic. Bulk packed reservoirs come in packs of 5 reservoirs.

Thomas No.	Description	Case
1228K15	Reagent Reservoir, 25 mL Divided, Sterile, Individually Wrapped	100
1228K16	Reagent Reservoir, 25 mL Divided, Sterile, Bulk Packed	200
1228K17	Reagent Reservoir, 25 mL Divided, Non-Sterile, Bulk Packed	100
1228K18	Reagent Reservoir, 25 mL, Sterile, Individually Wrapped	100
1228K19	Reagent Reservoir, 25 mL, Sterile, Bulk Packed	200
1228K20	Reagent Reservoir, 25 mL, Non-Sterile, Bulk Packed	100
1228K21	Reagent Reservoir, 50 mL, Sterile, Individually Wrapped	100
1228K22	Reagent Reservoir, 50 mL, Sterile, Bulk Packed	200
1228K23	Reagent Reservoir, 50 mL, Non-Sterile, Bulk Packed	100
1228K25	Reagent Reservoir, 100 mL, Sterile, Bulk Packed	200
1228K26	Reagent Reservoir, 100 mL, Non-Sterile, Bulk Packed	100