



RC 5

Heating Plates

Single-position heating plate with aluminum top



Excellent Heating Performance

The RC 5 is equipped with a 135 mm diameter aluminum alloy top plate ensuring a great temperature homogeneity and optimum heat transfer across the entire surface.

Designed to guarantee ease of use, the longest lifespan and excellent resistance to chemicals.

- Maximum temperature 310 °C
- Accepts flasks of up to 15 L volume
- Hot plate warning indicators
- IP 42 certified



INSTRUMENT - CODE

RC 5	230 V / 50-60 Hz	F20700610
RC 5	115 V / 50-60 Hz	F20710610

TECHNICAL DATA

	RC 5
HOUSING MATERIAL	Epoxy painted aluminum structure
PLATE MATERIAL	Aluminum alloy
PLATE DIMENSION	135 mm
TEMPERATURE RANGE	Room temp. - 310 °C
HOT TOP WARNING SYSTEM	Illuminated Icon (when the temperature is above 50 °C)
POWER INPUT	630 W
PERMISSIBLE AMBIENT TEMPERATURE	5 - 40 °C
PROTECTION CLASS	IP 42
WEIGHT	1.4 Kg (3.1 lb)
DIMENSIONS (WxHxD)	160x270x85 mm (6,3x10,6x3,3)



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
40, Burt Drive, Unit #1, Deer Park
NY 11729 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

LATAM
velplatam@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +8621 34500630
velpchina@velp.com