

## PDCD7 Polyclonal Antibody

Catalog number: 12485-1-AP


Size: 38 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PDCD7; ES18, HES18, PDCD7,  
programmed cell death 7



HeLa cells were subjected to SDS  
PAGE followed by western blot with  
12485-1-AP( PDCD7 Antibody) at  
dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PDCD7 MW:	485aa,55kd
Observed PDCD7 MW:	55 kDa
Positive WB detected in	HeLa cells, Jurkat cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3175
GenBank accession number:	BC016992
Gene ID (NCBI):	10081
Full name:	Programmed cell death 7

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.