

15111

Procedure Mask

- Complies with ASTM F2100 Level 1 requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning




Thomas No.	Mfr. No.	Description
22A00C289	15111	Procedure mask, ASTM Level 1, Earloop, Blue, BFE >99.6

Specifications

Product	Procedure Face Mask
Color	Blue
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.6
Filtration - PFE	≥98.7
Delta P	≤2.4
Flame Spread	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Case

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1